

Applause™ 3'-Amp System (5100-24)

Enter the following information to automatically calculate the volumes needed to prepare each reaction. The calculated volumes include an appropriate overfill in excess of the nominal volume requirements to allow for volume loss due to handling. Simply print this document to create a working guide for your experiment, which can be kept as a record.

Operator's Name: _____ Date: _____

Applause 3'-Amp System Lot No.: _____

SPIA® Product Purification Kit Name and Lot No.: _____ Number of Samples:* _____

Thermal Cycler Programs	
FIRST STRAND cDNA SYNTHESIS	
Program 1: Primer Annealing	65°C – 5 min, hold at 4°C
Program 2: First Strand Synthesis	4°C – 1 min, 48°C – 60 min, 70°C – 15 min, hold at 4°C
SECOND STRAND cDNA SYNTHESIS	
Program 3: Second Strand Synthesis	4°C – 1 min, 25°C – 10 min, 50°C – 30 min, 70°C – 5 min, hold at 4°C
SPIA AMPLIFICATION	
Program 4: SPIA® Amplification	4°C – 1 min, 47°C – 90 min, 95°C – 5 min, hold at 4°C

* Number of samples field ties into embedded logic to calculate master mix volumes
NuGEN recommends processing a minimum of 8 samples at a time.

First Strand cDNA Synthesis			
Thaw the First Strand Reagents (blue) and Nuclease Free Water D1 (green) . Mix each reagent, spin and place on ice.			
Add 2 μ L of A1 into each 0.2 mL PCR tube.			
For each sample place 5 μ L of RNA (50–200 ng) into a 0.2 mL PCR tube, mix and spin.			
Place the tubes in a thermal cycler running Program 1 (65°C – 5 min, hold at 4°C).			
Prepare First Strand Master Mix (calculation allows for appropriate overfill). Please be sure to pipet A3 enzyme slowly and rinse out tip at least five times into buffer. Per sample combine: 2.5 μ L Buffer Mix A2 VER 6 + 0.5 μ L Enzyme Mix A3 VER 1 Mix well.	No. of Samples	A2	A3
	1	2.5 μ L	0.5 μ L
Add 3 μ L of the First Strand Master Mix to each tube, mix and spin.			
Place the tubes in a thermal cycler running Program 2. (4°C – 1 min, 48°C – 60 min, 70°C – 15 min, hold at 4°C)			

Second Strand cDNA Synthesis			
Thaw the Second Strand Reagents (yellow) . Vortex B1 . Flick B2 to mix. Spin all and place on ice.			
Once the thermal cycler reaches 4°C, remove tubes, spin and place on ice.			
Prepare Second Strand Master Mix (calculation allows for appropriate overfill). Please be sure to pipet B2 enzyme slowly and rinse out tip at least five times into buffer. Per sample combine: 9.5 μ L Buffer Mix B1 VER 3 + 0.5 μ L Enzyme Mix B2 VER 2. Mix well.	No. of Samples	B1	B2
	1	9.5 μ L	0.5 μ L
Add 10 μ L of Second Strand Master Mix to each first strand reaction tube, mix and spin.			
Place the tubes in a thermal cycler running Program 3 (4°C – 1 min, 25°C – 10 min, 50°C – 30 min, 70°C – 5 min, hold at 4°C).			
Once the thermal cycler reaches 4°C, spin and place tubes on ice.			

SPIA Amplification				
Thaw the SPIA Amplification Reagents (red) . Vortex C2 and C1 , invert C3 5 times. Spin all, place on ice.				
Prepare SPIA Master Mix (calculation allows for appropriate overfill). Per sample combine: 36 μ L C2 VER 6 + 4 μ L C1 VER 10 + 20 μ L C3 VER 5 . Mix well.	No. of Samples	C2	C1	C3
	1	36 μ L	4 μ L	20 μ L
Add 60 μ L of the SPIA Master Mix to each Second Strand synthesis reaction tube, mix thoroughly and spin down.				
Place tubes in a thermal cycler running Program 4 (4°C – 1 min, 47°C – 90 min, 95°C – 5 min, hold at 4°C).				
Once the thermal cycler reaches 4°C, spin and place tubes on ice.				
Proceed immediately to purification step or store SPIA cDNA at -20°C.				

Purification of Amplified SPIA cDNA		
Refer to the user guide and follow the method of choice for purification:	Purification Kit Part No.	Purification Kit Lot No.
Add Binding Buffer in volume of:	Spin at speed:	For a duration of:
Add Wash Buffer in volume of:	Spin at speed:	For a duration of:
Repeat for second wash.		
To elute sample use Nuclease-Free Water D1 provided with the Ovation Kit.		
Add Nuclease-Free Water D1 in volume of:	Spin at speed:	For a duration of: